capturer

EN	INSTRUCTIONS FOR USE THROMBUS EXTRACTION CATHETER	2
ES	INSTRUCCIONES DE USO. CATÉTER EXTRACTOR DE TROMBOS	7
PT	INSTRUÇÕES DE UTILIZAÇÃO CATETER EXTRATOR DE TROMBOS	12
IT	ISTRUZIONI PER L' USO CATETERE PER ESTRAZIONE DEI TROMBI	17
FR	INSTRUCTIONS D'UTILISATION CATHÉTER D'EXTRACTION DE THROMBUS	22
DE	GEBRAUCHSANWEISUNG THROMBUS-EXTRAKTIONSKATHETER	27
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ КАТЕТЕР ДЛЯ АСПИРАЦИИ ТРОМБОВ	32
CS	NÁVOD K POUŽÍTÍ EXTRAKČNÍ KATETR TROMBU	37
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ Καθετήρας αναρρόφησης θρόμβου	42
PL	CEWNIK DO ODSYSANIA SKRZEPLIN INSTRUKCJA OBSŁUGI	47
SE	INSTRUKTIONER FÖR ANVÄNDNING AV KATETER FÖR EXTRAKTION AV TROMBOSER	52
TU	TROMBUS ÇIKARMA KATETERİ KULLANIM TALİMATLARI	57
NL	INSTRUCTIES VOOR TROMBUS EXTRACTIEKATHETER.....	62


INSTRUCTIONS FOR USE THROMBUS EXTRACTION CATHETER

1. Device description

The Capturer thrombus extraction catheter is a rapid exchange catheter, having a single proximal lumen and a double distal lumen.

The catheter is compatible with 0.014 inch /0.36 mm guidewires.

The larger lumen runs the length of the catheter from its proximal end to its distal part and it is intended for the removal of the thrombus when connected to one of the syringes included in the system, while the smaller lumen is only in the last distal 17,5 cm of the catheter, and it is able to accommodate the guidewires facilitating advancement of the catheter to the area that must be treated.

The proximal region is stiff to allow for the catheter to be pushed, whereas the distal region is flexible to allow for it to navigate within the artery. The catheter's distal region is coated by a durable hydrophilic coating which confers lubricity for the catheter to navigate through the arteries.

A radiopaque marker is placed inside the tip of the catheter to indicate, at all times, the position of the catheter within the artery. At 15 mm of the radiopaque marker, there is another radiopaque marker to provide an idea of how long the injury is. In the distal region before the tip, there is a hole through which the thrombus can be aspirated by connecting the device to the included manual extraction syringe.

The extraction catheter is equipped with a stiffening mandrel (stylet) coupled to the catheter's luer connector. The use of this mandrel (stylet) is optional, if you do not wish to introduce the catheter with stiffening mandrel (stylet), it is removed before the procedure begins. The purpose of this stiffening mandrel (stylet) is to prevent the catheter from kinking during its advancement towards the treatment area. Before beginning aspiration, the mandrel (stylet) must be removed to maximize all the aspiration area.

Two extraction syringes, a stopcock, an extension line and two filter baskets to collect and deposit the clots removed from the patient are included.

The catheter has an useful length of 140 cm, it is provided under two reference numbers:

- DET R14 145 601 - 6F Thrombus extraction catheter: compatible to be introduced inside a 6F guide catheter with minimum internal diameter of 1.78 mm.
- DET R14 145 701 - 7F Thrombus extraction catheter: compatible to be introduced inside a 7F guide catheter with minimum internal diameter of 2.00 mm.

Contents

- One iVascular Capturer thrombus extraction catheter with an optional stiffening mandrel
- Two filter baskets
- One stopcock
- One extension line
- Two vacuum syringes with locking plunger
- One leaflet with instructions for use

2. Indications

The Capturer thrombus extraction catheter is especially indicated for the removal of fresh, soft emboli and thrombi from the coronary or peripheral vasculature, with a minimum diameter of 2 mm.

3. Contraindications

- Vessels less than 2mm in diameter.

INSTRUCTIONS FOR USE THROMBUS EXTRACTION CATHETER

- Coronary artery spasm without significant stenosis.
- Extraction of fibrous, adherent or calcified material, such as chronic clots or atherosclerotic plaques.
- The venous system.
- Contraindications to anticoagulant or antiplatelet treatment.

4. Warnings

- The contents are supplied STERILE. Inspect the packaging before opening it. If any defects are found or if it is damaged, discard the product and return it to the manufacturer or distributor.
- Check the expiration date and do not use products that have surpassed this date.
- The product is meant to be used a single time. It must not be re-sterilized or reused once the procedure has been concluded. Reutilization of the device in another patient may cause cross contamination, infections or transmission of infectious diseases from one patient to the next. Reutilizing the device may cause it to alter and limit its effectiveness.
- Use aseptic techniques when the device is removed from its primary packaging. Do not wipe dry with gauze.
- Do not expose the device to organic solvents, or contrast media not indicated for its intravascular use.
- Once inside the vascular system, the catheter must be handled under high-quality fluoroscopic vision. If any resistance is met, determine cause before continuing with procedure or remove catheter.
- Advance the device on the guidewire with fluoroscopic vision. Do not permit the advancement of the device without the guidewire inside it.
- Do not use the catheter for the delivery or infusion of diagnostic, embolic or therapeutic materials into blood vessels, as it has been designed for other uses.
- Do not perform high pressure contrast injections through the Capturer thrombus aspiration catheter while using a guide catheter of the same French size, as it could damage it, making it difficult to remove from the guide catheter of the same French size.
- Once the thrombus extraction has started through the extraction lumen, do not attempt to infuse solutions as possible thrombi that may have remained in the catheter, may be released again into the vasculature, which may cause thromboembolic events, serious injury and even death. Remove the catheter, clean it and introduce it again if necessary.
- If during thrombus extraction, flow into the syringe stops or is restricted, do not attempt to flush the extraction lumen while it is inside the patient. Remove the catheter from the patient, flush the extraction lumen cleaning it carefully and reinsert it into the patient if necessary.
- Percutaneous intervention must be performed in hospitals equipped with emergency media to perform open-heart operations or in hospitals with an easy access to other hospitals that can perform this surgery.

5. Precautions

- The device should be used by experienced physicians thoroughly trained who know well the percutaneous, transluminal angioplasty technique.
- Administer the adequate medical therapy to the patient: anticoagulants, vasodilators, etc, in accordance with the insertion protocol of intravascular catheters.
- Inspect the catheter prior to use for any bends or kinks and do not use any damaged catheter.
- Excessive tightening of the hemostatic valve onto the catheter may result in damage.
- Exercise the utmost care in handling the catheter so as not to damage it during its advancement.
- In case the stiffening mandrel (stylet) is used, remove it before beginning aspiration otherwise, the aspiration area is decreased.
- When using the extraction catheter with a guide catheter of the same French size, there will be insufficient room to deliver contrast fluid around the guide catheter.
- Check that all the fittings are secured firmly to prevent air from being introduced into the system during thrombus extraction, as it would decrease effectiveness of extraction.

INSTRUCTIONS FOR USE

THROMBUS EXTRACTION CATHETER

- Store in a cool, dry place, away from direct sunlight.
- After its use, this device may pose a biological hazard. Handle it and dispose of it in compliance with accepted medical practices and the pertinent local, state or federal laws and regulations.

6. Possible adverse effects/complications

Among the possible adverse effects and/or complications that might appear before, during or after procedure, the following may be cited:

- Death
- Acute myocardial infarction
- Unstable Angina
- Ictus
- Distal clot and /or plaque embolization
- Thrombosis
- Perforation and rupture of blood vessels
- Artery spasms
- Arrhythmias
- Total arterial occlusion
- Hematomas
- Hemorrhage
- Infections
- Arteriovenous Fistula
- Changes in blood pressure
- Sensitivity or allergic reaction to contrast fluid, drugs and materials

7. Instructions for use

7.1. Necessary components

- Normal heparinized saline
- Contrast medium (use contrast media indicated for intravascular use)
- 0.014 guidewire
- 6F/7F Introducer
- 6F/7F Guiding catheter
- Stopcock
- Extraction syringes
- Filter basket
- Several 10-20 cc standard syringes with saline to flush the system

7.2. Preparation of the catheter

Before beginning the procedure, carefully examine all the equipment to be used. Look out for possible damage in product packaging and discard those that are not in perfect conditions. These instructions only provide technical assistance, but do not obviate the need for formal training.

- Open the headerbag and empty the tray. Transfer the contents into the sterile field. Remove the stylet and inspect the catheter to make sure it is in good condition.
- Steep the catheter in sterile saline to activate hydrophilic coating and improve the catheter's ability to navigate.
- Connect the extension line to the luer connection of the catheter and after that the stopcock. Leave the stopcock in the "OFF" position.

NOTE: In case the stiffening mandrel (stylet) is used, connect the extension to the mandrel's luer lock.

- Flush the inside of the catheter by attaching a syringe with 10 ml of sterile saline to the extension line and injecting it through the system (to do so turn on the stopcock and off again when the process has finished), checking that the fluid comes out of the distal hole and that there are no

INSTRUCTIONS FOR USE THROMBUS EXTRACTION CATHETER

leaks. Remove syringe.

7.3. Insertion technique/treatment

- Insert the introducer and the guiding catheter by means of standard techniques recommended in the manufacturer's instructions for use.
- Insert the angioplasty guidewire through the guiding catheter and place it in the injury using standard techniques recommended in the manufacturer's instructions for use.
- Before introducing the guidewire in the distal end of the catheter, flush the lumen with sterile saline. Immediately after introduce the guidewire from the distal end making sure that the guide exits through the exit port located approximately 15 cm of the distal end. In order to perform this operation it may be necessary to straighten the catheter to facilitate the guidewire's passage. Stretch the guidewire proximally keeping its distal end within the catheter for its protection.
- Open the hemostatic valve placed on the guiding catheter connector, and advance the extraction catheter approximately 40 cm before closing the valve adaptor. Prevent the hemostatic valve adaptor from damaging the catheter shaft which may influence later on the catheter's properties, such as its extraction capacity.
- Advance the catheter over the guidewire. If any resistance is found, stop advancement of the catheter and do not continue until the cause of the resistance has been determined by fluoroscopy.
- Continue to advance the extraction catheter over the guidewire to the tip of the guiding catheter. Try not to rotate or bend the catheter around the guiding catheter to avoid damaging it.
- Once outside the guiding catheter, the extraction catheter may be easily observed thanks to the distal radiopaque mark, which indicates the catheter's advancement and position. Advance the catheter through the angioplasty guide to the thrombus area.
- Once the thrombus area is reached, remove the stiffening mandrel (stylet) and connect the extension directly to the aspiration catheter's luerlock connector. Dip the stiffening mandrel (stylet) in heparinized saline to clean it. It might be necessary to reintroduce it into the catheter in case the procedure is repeated.
- Fill the extraction syringe with 2 ml of heparinized saline. Connect the syringe through the stopcock (which must be in the "OFF" position).

In the extraction, proceed as follows:

- With the stopcock in the "OFF" position (closed), pull back the plunger until the desired extraction volume mark has been reached. Twist the plunger to lock the syringe in the vacuum position. No bubbles must pass through the syringe, otherwise, it might mean that there is a leak in the connection between the stopcock and the syringe.
- Check the catheter's position fluoroscopically and open the stopcock to begin extraction. Advance the catheter distally. Blood will enter the syringe until all of the vacuum is gone. Should aspiration not fill the syringe, remove the catheter from the patient, flush the extraction lumen to clean it thoroughly and place it again in vasculature.
- Once the extraction procedure is finished, turn the stopcock in the "OFF" position again and begin removing the catheter or connect the other syringe to continue extraction.
- Remove the catheter completely. Keep the guidewire in its position.
- Perform an angiography through the guiding catheter to estimate if it is necessary to extract more thrombi. If it is, flush the inside of the extractor making sure that no residues are left or use a new catheter and perform the entire procedure again just as described in the previous paragraphs.
- Once the procedure is finished, remove the thrombus extracting device. Maintain the guiding catheter and the guidewire in case of proceeding with a stent or other device implantation.
- When all the procedures have concluded and the devices have been removed, leave the introducer in situ until the hemodynamic profile returns to normal.
- Blood and thrombi extracted into the syringe may be filtered for subsequent laboratory analysis.

INSTRUCTIONS FOR USE THROMBUS EXTRACTION CATHETER

using the filter basket in the package.

Warranty

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

STANDARDIZED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES (EN 980:2008 STANDARD)

Catalogue reference

Batch no.

Expiry date

Keep away from sunlight

Manufacturer

Sterilized with ethylene
oxide

Limitation in storage
temperature

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