oceanus

18

EN	INSTRUCTIONS FOR USE. BALLOON CATHETER FOR PERIPHERAL PREDILATION, FOR 0.018" GUIDE WIRE, WITH HYDROPHILIC COATING	3
ES	INSTRUCCIONES DE USO. CATÉTER BALÓN DE PREDILATACIÓN PERIFÉRICO, PARA GUÍA DE ALAMBRE DE 0.018", CON RECUBRIMIENTO HIDROFÍLICO	7
PT	INSTRUÇÕES DE USO. CATETER BALÃO DE PRÉ-DILATAÇÃO PERIFÉRICO, PARA FIO-GUIA DE 0,018", COM REVESTIMENTO HIDROFÍLICO	11
IT	ISTRUZIONI PER L'USO. CATETERE A PALLONCINO PER DILATAZIONE PERIFERICA (PTA), CON RIVESTIMENTO IDROFILICO E COMPATIBILE CON FILO GUIDA DA 0.018"	15
FR	MODE D'EMPLOI. CATHÉTER A BALLONNET DE PRÉ-DILATATION PÉRIPHÉRIQUE, POUR FIL GUIDE DE 0,018", AVEC REVÊTEMENT HYDROPHILE	19
DE	GEBRAUCHSHINWEISE. BALLONKATHETER FÜR PERIPHERE VORDILATATION, FÜR 0,018" FÜHRUNGSDRAHT MIT HYDROPHILER BESCHICHTUNG	23
RU	ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ. БАЛЛОННЫЙ КАТЕТЕР С ГИДРОФИЛЬНЫМ ПОКРЫТИЕМ ДЛЯ ПРЕДИЛАТАЦИИ ПЕРИФЕРИЧЕСКИХ СОСУДОВ, ДЛЯ ПРОВОДНИКА 0,018"	27
CS	NÁVOD K POUŽITÍ. BALÓNKOVÝ KATÉTR PRO PERIFERNÍ PREDILATACI S VODICÍM DRÁTEM 0,018", S HYDROFILNÍM POUŠKEM	31
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ. ΚΑΘΗΤΡΑΣ ΜΠΑΛΟΝΙ ΓΙΑ ΠΡΟΔΙΑΣΤΟΛΗ ΠΕΡΙΦΕΡΙΚΩΝ ΑΡΤΗΡΙΩΝ ΓΙΑ ΣΥΡΜΑΤΙΝΟ ΟΔΗΓΟ 0,018". ΜΕ ΥΔΡΟΦΙΛΗ ΕΠΙΚΑΛΥΨΗ	35
PL	INSTRUKCJA UŻYTKOWANIA. CEWNIK BALONOWY DO PREDYLATACJI NACZYŃ OBWODOWYCH, DO PROWADNIKÓW 0,018", Z POWŁOKĄ HYDROFILINĄ	39

INSTRUCTIONS FOR USE
BALLOON CATHETER FOR PERIPHERAL PREDILATION,
FOR 0.018" GUIDE WIRE, WITH HYDROPHILIC COATING

1. Description of the Product

The balloon catheter is a coaxial catheter (dual lumen from the connector to the tip, also known as OTW), designed for percutaneous transluminal angioplasty of peripheral arteries.

The catheter has a coaxial dual-lumen body from the connector to the tip.

On the distal part of the catheter, before the tip, is located the balloon or inflatable segment, to different pressures. The balloon dilates the artery when is inflated by means of infusion of contrast medium inside it. The connector is Y-shaped, with two ports:

- The side port allows contrast medium passage to dilate The balloon.
- The straight port is intended for guide wire insertion. The guide wire is inserted through the tip of the catheter and comes out by the connector port to lead the catheter to the lesion. The guidewire may also be inserted from the coaxial luer connector and go through the catheter to the tip.

Moreover, the lumen inner of the catheter is the one through which the guide wire runs, while the outer lumen is intended for contrast medium passage.

The maximum diameter of the guide wire must not exceed 0.46mm = 0.018". Alternatively, it can be used a guide wire of 0.014" = 0,36mm.

The balloon is intended to reach different diameters at different pressures (according to the compliance curve included on the primary packaging).

To inflate the balloon, the connector inflation port must be connected to an inflation device. The balloon is inflated to predictable diameters with the pressure shown by the manometer. Two radiopaque markers are located at each end of the balloon in order to mark its length and help the user to see the catheter while inside the patient.

At the distal end is the tip, which is rounded and atraumatic in shape in order to avoid damaging the arteries while it is being advanced.

The distal tube of the catheter is coated with a durable hydrophilic coating to help the catheter go through tortuous arteries smoothly.

The useful length of the catheter is 100 cm, 140 cm or 150 cm depending on the item reference.

The balloon catheter is supplied in the following balloon diameters and lengths:

BALLOON DIMENSIONS		Lengths (mm)						
		20	40	60	80	120	150	200
Diameters (mm)	2,0	X	X	X	X	X	X	X
	2,5	X	X	X	X	X	X	X
	3,0	X	X	X	X	X	X	X
	3,5	X	X	X	X	X	X	X
	4,0	X	X	X	X	X	X	X
	5,0	X	X	X	X	X	X	X
	6,0	X	X	X	X	X	X	X
	7,0	X	X	X	X	X		
8,0	X	X	X	X				

Contents

- One OTW catheter covered by protection sheath for the balloon and a stylet to protect the guide wire lumen. This is inserted into a dispenser in a sterile bag.
- One card with the compliance curve showing the nominal inflation pressure and the recommended maximum pressure.
- One leaflet with instructions for use.

INSTRUCTIONS FOR USE
BALLOON CATHETER FOR PERIPHERAL PREDILATION,
FOR 0.018" GUIDE WIRE, WITH HYDROPHILIC COATING

2. Indications

The device is indicated to dilate stenoses in iliac, femoral, popliteal, infra-popliteal and renal arteries and for treating obstructive lesions of native or synthetic arteriovenous dialysis fistulae, with a reference diameter balloon from 2 and 8 mm and lengths from 20 up to 200 mm.

3. Contraindications

Inability to cross lesion with guide-wire.

4. Warnings

- The device should only be used by experienced doctors familiar in the percutaneous transluminal angioplasty (PTA) technique.
- The percutaneous procedure must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospitals where this type of operation can be performed.
- The product is for single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the packaging before opening it. If there are any defects or the packaging is damaged, do not use the product.
- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter according to the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than the artery to be treated.
- Do not dry with gauze.
- Do not expose the product to organic solvents or to contrast media not suitable for intravascular use.
- Do not use air or gaseous media to inflate the balloon.
- Do not continue with the procedure if there is friction.
- Advance the product on the guide wire using fluoroscopy. Do not allow the product to advance without the guide wire inside it.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.

5. Precautions

- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the procedure for inserting intravascular catheters.
- Use with standard 0.018" guide wires (alternatively 0.014 ") and standard introducers. If the introducer length goes above 25cm or the braided type is used, it may be necessary to increase its size to avoid catheter friction.
- Before use, check the size and suitability of the catheter for its intended use.
- Proceed with extreme care so as not to damage the catheter while advancing it. If during advance there is any friction stop and determine the cause before continuing.
- The balloon must be inflated with a mixture of saline solution and contrast medium (preferably 50/50, though the proportion of contrast medium can be reduced in large balloons, like 4mm diameter or longer ones).
- If any resistance is noticed during advance, stop and determine its cause before continuing.
- Do not exceed the rated burst pressure (RBP) as the balloon could burst.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guide wire and introducer be removed as a single unit.
- Store in a cool, dry place away from direct sunlight.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.

INSTRUCTIONS FOR USE
BALLOON CATHETER FOR PERIPHERAL PREDILATION,
FOR 0.018" GUIDE WIRE, WITH HYDROPHILIC COATING

6. Possible Adverse Effects / Complications

The possible adverse effects and/or complications which might arise before, during or after the procedure include the following:

- Death
- Ictus / embolism / thrombosis
- Haemodynamic deterioration
- Acute elastic recoil
- Arrhythmia
- Aneurysm or pseudo-aneurysm
- Infections
- Total occlusion of the artery
- Spasm
- Perforation or dissection of the treated area
- Arteriovenous fistula
- Local bleeding with haematoma at the access site
- Hypo/hypertension
- Allergic reactions to materials.

7. How to Use

7.1. Equipment Required

- Heparinized normal saline solution.
- Contrast medium (use contrast media suitable for intravascular use).
- 0.018" guidewire. Alternatively 0.014".
- Introducer with a haemostatic valve of The size specified on The label. Do not use smaller sizes as these might damage The catheter.
- Three-way stopcock.
- Line extensions.
- Inflation device (manual pump with built-in manometer).
- Several standard 10-20cc syringes with saline solution for washing the system.

7.2. Preparing the Catheter

- Remove the catheter from the protective dispenser. Examine the catheter, if there is any anomaly do not use it. Check that it is the right size. Gently remove the protective sheath and stylet from the balloon."Connect a 10ml syringe containing sterile saline solution to the port for the guide wire (straight part of the connector) and irrigate the lumen until the liquid comes out of the tip.
- Purge the air from the whole system, connectors and pump before starting inflation:
 - Attach a three-way stopcock to the contrast liquid port (the side part of the connector). An extension line could be optionally placed between the balloon connector and the three-way stopcock.
 - Close the air passage through the balloon.
 - Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
 - Open the passage in the three-way stopcock between the syringe and the catheter.
 - With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
 - When bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe.

DANGER: If bubbles do not stop coming from the catheter into the syringe after 3 minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly closed. If bubbles continue

INSTRUCTIONS FOR USE
BALLOON CATHETER FOR PERIPHERAL PREDILATION,
FOR 0.018" GUIDE WIRE, WITH HYDROPHILIC COATING

to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

- In order to activate the hydrophilic coating, it is recommended to wet the catheter with a sterile saline solution immediately prior to its insertion through the introducer sheath.

7.3. Technique for Insertion / Treatment

- Position the 0.018" guide wire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times. If previously placed a 0.014" guide wire there is no need to remove it, as the catheter is also compatible with that size.
- Insert the proximal end of the guide wire into the distal tip of the balloon catheter.
- Carefully advance the balloon catheter on the guide wire, through the introducer and with the aid of fluoroscopy, until it reaches the area to be treated. Do not allow advance of the product without guidewire inside. In the event that the balloon catheter selected cannot be placed inside the lesion, take a smaller size to predilate the lesion.
- Connect the inflation device (1/3 of its capacity of contrast liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Apply pressure for 30-60 seconds depending on the lesion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 and 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire in position, carefully retract the catheter until it can be removed through the introducer.
- Inject contrast liquid through the introducer to check the treated lesion before removing the guide wire.

8. Warranty

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

STANDARDIZED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES
(EN 980:2008 STANDARD)

			
Catalogue reference	Batch no.	Expiry date	Keep away from sunlight
			
Manufacturer	Sterilized with ethylene oxide	Limitation in storage temperature	Keep dry
			
Consult instructions for use	Do no re-use	Do not re-sterilize	Do not use if the pack is damaged

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